


Anti-Pomp Rabbit Polyclonal Antibody For Wb, Elisa - Q9Y244

Item Number: CM0029773



Introduction

Rabbit polyclonal antibody against human POMP, validated for Western blot and ELISA. Ideal for proteasome maturation studies, chaperone research, and protein degradation analysis.

[Learn More](#)

Parameter	Value
Product Name	POMP Rabbit pAb
Remarks/Alias	UMP1; PRAAS2; HSPC014; C13orf12; PNAS-110; POMP
Species	Human
GeneID (Human)	51371
GeneID	51371
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-141 of human POMP (NP_057016.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q9Y244
Protein Weight	16kDa
Shipping	Ice bag