

## Anti-Pou2F1/Oct1 Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - P14859

Item Number: CM0008432



### Introduction

Rabbit monoclonal antibody targeting human POU2F1/OCT1, validated in WB, IHC-P, and ELISA. Detects human, mouse, and rat POU2F1 (~76 kDa). Suitable for studies on transcription regulation, nuclear receptor signaling, and HSV infection mechanisms.

[Learn More](#)

| Parameter        | Value                                  |
|------------------|----------------------------------------|
| Product Name     | POU2F1/OCT1 Rabbit mAb                 |
| Remarks/Alias    | OCT1; OTF1; Oct1Z; oct-1B; POU2F1/OCT1 |
| Species          | Human                                  |
| GeneID (human)   | 5451                                   |
| GeneID           | 5451                                   |
| Immunogen        | Recombinant protein                    |
| Source           | Rabbit                                 |
| Category         | Monoclonal Antibodies                  |
| Application      | WB, IHC-P, ELISA                       |
| Cross Reactivity | Human, Mouse, Rat                      |
| SWISS            | P14859                                 |
| Protein Weight   | 76 kDa                                 |
| Shipping         | Ice bag                                |