

Anti-Pou2F1/Oct1 Rabbit Polyclonal Antibody For Wb, Elisa - P14859

Item Number: CM0022262



Introduction

Rabbit polyclonal antibody against human POU2F1/OCT1, offered for WB and ELISA applications, cross-reacts with rat, with a predicted molecular weight of 76 kDa. Ideal for transcription factor and chromatin studies. UniProt P14859.

[Learn More](#)

Parameter	Value
Product Name	POU2F1/OCT1 Rabbit pAb
Remarks/Alias	OCT1; OTF1; Oct1Z; oct-1B; POU2F1/OCT1
Species	Human
GeneID (Human)	5451
GeneID	5451
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 15-290 of human POU2F1/OCT1 (NP_001185712.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Rat
SWISS	P14859
Protein Weight	76kDa
Shipping	Ice bag