

Anti-Pou3F1 Polyclonal Antibody For Wb, Elisa - Q03052

Item Number: CM0022417

Introduction

POU3F1 rabbit polyclonal antibody targeting human OCT6/SCIP transcription factor. Validated for Western blot and ELISA applications with cross-reactivity to mouse and rat. Ideal for neuroscience and stem cell research.

[Learn More](#)

Parameter	Value
Product Name	POU3F1 Rabbit pAb
Remarks/Alias	OCT6; OTF6; SCIP; POU3F1
Species	Human
Gene ID (Human)	5453
Gene ID	5453
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-260 of human POU3F1 (NP_002690.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q03052
Protein Weight	45kDa
Shipping	Ice bag