


Anti-Ppfia1 Polyclonal Antibody For Wb, Elisa - Q13136

Item Number: CM0013071



Introduction

Anti-PPFIA1 rabbit polyclonal antibody for Western blotting (WB) and ELISA. Validated in human samples, detects liprin-alpha-1 (UniProt Q13136, 136kDa) for cell adhesion and signaling research.

[Learn More](#)

Parameter	Description
Product Name	PPFIA1 Rabbit pAb
Remarks/Alias	LIP1; LIP.1; LIPRIN; PPFIA1
Species	Human
GeneID (Human)	8500
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 660-800 of human PPFIA1 (NP_003617.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q13136
Protein Weight	136kDa
Shipping	Ice bag