

Anti-Ppp4R3B Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q5Miz7

Item Number: CM0024526

Introduction

Rabbit polyclonal antibody targeting PPP4R3B (SMEK2), a regulatory subunit of PP4 localized to centrosomes/nucleus. Validated for WB, IF/ICC, ELISA with human, mouse, rat reactivity.

[Learn More](#)

Parameter	Value
Product Name	PPP4R3B Rabbit pAb
Remarks/Alias	PSY2; smk1; FLFL2; SMEK2; PP4R3B; PPP4R3B
Species	Human
GeneID (Human)	57223
GeneID	57223
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 628-817 of human PPP4R3B (NP_001269779.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q5MIZ7
Protein Weight	97kDa
Shipping	Ice bag