

Anti-Prkab2 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - O43741

Item Number: CM0025051



Introduction

High-quality rabbit polyclonal antibody targeting human PRKAB2 (AMPK beta-2 subunit). Validated for WB, IF/ICC, and ELISA applications with cross-reactivity to human, mouse, and rat. Ideal for metabolic and energy research.

[Learn More](#)

Parameter	Value
Product Name	PRKAB2 Rabbit pAb
Remarks/Alias	PRKAB2
Species	Human
GeneID (Human)	5565
GeneID	5565
Immunogen	Recombinant protein / Recombinant fusion protein containing a sequence corresponding to amino acids 1-272 of human PRKAB2 (NP_005390.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O43741
Protein Weight	30kDa
Shipping	Ice bag