


Anti-Prkcq (Npkc-Theta) Rabbit Polyclonal Antibody For Wb, Elisa - Q04759

Item Number: CM0030323



Introduction

Rabbit polyclonal antibody against human PRKCQ (nPKC-theta) for Western blot (WB) and ELISA. Cross-reacts with mouse and rat. Suitable for T-cell receptor signaling and immune response research.

[Learn More](#)

Parameter	Value
Product Name	PRKCQ Rabbit pAb
Remarks/Alias	PRKCT; nPKC-theta; PRKCQ
Species	Human
GeneID (Human)	5588
GeneID	5588
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-230 of human PRKCQ (NP_006248.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q04759
Protein Weight	82kDa
Shipping	Ice bag