

Anti-Prkcsh Polyclonal Antibody For Wb And Elisa - P14314

Item Number: CM0029663



Introduction

This rabbit polyclonal antibody reacts with human, mouse, and rat PRKCSH (Glucosidase 2 subunit beta) and is validated for Western blot and ELISA applications. Suitable for glycoprotein processing and polycystic kidney disease research.

[Learn More](#)

Parameter	Value
Product Name	PRKCSH Rabbit pAb
Remarks/Alias	GIIB; PCLD; PLD1; G19P1; PCLD1; PKCSH; AGE-R2; VASAP-60; PRKCSH
Species	Human
GeneID	5589
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 130-275 of human PRKCSH (P14314)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P14314
Protein Weight	59kDa
Shipping	Ice bag