


Anti-Prmt5 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - O14744

Item Number: CM0015493



Introduction

Rabbit polyclonal antibody against human PRMT5 (O14744), validated for WB, IHC-P, IF/ICC, ELISA. Cross-reacts with mouse and rat. Suitable for studies on arginine methylation, chromatin regulation, and spliceosome assembly.

[Learn More](#)

Parameter	Value
Product Name	PRMT5 Rabbit pAb
Remarks/Alias	HSL7; JBP1; SKB1; IBP72; SKB1Hs; HRMT1L5; PRMT5
Species	Human
GeneID (Human)	10419
GeneID	10419
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-250 of human PRMT5 (NP_006100.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O14744
Protein Weight	73kDa
Shipping	Ice bag