

Anti-Prolactin-Inducible Protein (Pip) Polyclonal Antibody For If/Icc, Elisa - P12273

Item Number: CM0027896



Introduction

Rabbit polyclonal antibody against human Prolactin-inducible protein (PIP/GCDFP-15). Validated for IF/ICC and ELISA. Recognizes 17 kDa mature secreted glycoprotein. Cross-reacts with human. For breast cancer and exocrine gland research.

[Learn More](#)

Parameter	Value
Product Name	PIP Rabbit pAb
Remarks/Alias	GPIP4; BRST-2; GCDFP15; GCDFP-15; PIP
Species	Human
Gene ID (Human)	5304
Gene ID	5304
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 29-146 of human PIP (NP_002643.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	IF/ICC, ELISA
Cross Reactivity	Human
SWISS	P12273
Protein Weight	17kDa
Shipping	Ice bag