


## Anti-Prp31 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q8Wwy3

Item Number: CM0029294



### Introduction

Validate pre-mRNA splicing with this PRPF31 rabbit polyclonal antibody. Specific for human, mouse, rat; ideal for WB, IHC-P, ELISA. 55 kDa U4/U6 snRNP component essential for spliceosome assembly.

[Learn More](#)

| Parameter        | Value                                                                                                              |
|------------------|--------------------------------------------------------------------------------------------------------------------|
| Product Name     | PRPF31 Rabbit pAb                                                                                                  |
| Remarks/Alias    | RP11; PRP31; SNRNP61; NY-BR-99; PRPF31                                                                             |
| Species          | Human                                                                                                              |
| Gene ID (Human)  | 26121                                                                                                              |
| Gene ID          | 26121                                                                                                              |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 1-200 of human PRPF31 (NP_056444.3). |
| Source           | Rabbit                                                                                                             |
| Category         | Polyclonal Antibodies                                                                                              |
| Application      | WB, IHC-P, ELISA                                                                                                   |
| Cross Reactivity | Human, Mouse, Rat                                                                                                  |
| SWISS            | Q8WWY3                                                                                                             |
| Protein Weight   | 55kDa                                                                                                              |
| Shipping         | Ice bag                                                                                                            |