

Anti-Prr12 Rabbit Polyclonal Antibody For Wb, Elisa - Q9ULL5

Item Number: CM0020476

Introduction

Rabbit polyclonal antibody targeting human PRR12, validated for western blot and ELISA with cross-reactivity to mouse and rat. Suitable for chromatin biology and neurodevelopment research. - Q9ULL5

[Learn More](#)

Parameter	Value
Product Name	PRR12 Rabbit pAb
Remarks/Alias	NOC; KIAA1205; PRR12
Species	Human
Gene ID (Human)	57479
Gene ID	57479
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1687-2036 of human PRR12 (NP_065770.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9ULL5
Protein Weight	211kDa
Shipping	Ice bag