

Anti-Prss1 Polyclonal Antibody For Wb, Ihc-P, Elisa - P07477

Item Number: CM0015965



Introduction

Rabbit polyclonal antibody targeting human PRSS1 (trypsin-1). Validated for WB, IHC-P, and ELISA. Recognizes human and rat PRSS1. Suitable for research on serine protease digestion and extracellular proteolysis.

[Learn More](#)

Parameter	Details
Product Name	PRSS1 Rabbit pAb
Remarks/Alias	TRP1; TRY1; TRY4; TRYP1; PRSS1
Species	Human
GeneID	5644
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 24-247 of human PRSS1 (NP_002760.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Rat
SWISS	P07477
Protein Weight	24kDa
Shipping	Ice bag