

## Anti-Prss2 Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - P07478

Item Number: CM0002412



### Introduction

Anti-PRSS2 rabbit monoclonal antibody validated for WB, IHC-P, and ELISA. Recognizes human, mouse, and rat trypsin-2. Ideal for studying digestive proteases and defensin processing in Paneth cells.

[Learn More](#)

| Parameter        | Value                                                                                                               |
|------------------|---------------------------------------------------------------------------------------------------------------------|
| Product Name     | PRSS2 Rabbit mAb                                                                                                    |
| Remarks/Alias    | TRY2; TRY8; TRYP2; PRSS2                                                                                            |
| Species          | Human                                                                                                               |
| Gene ID (Human)  | 5645                                                                                                                |
| Gene ID          | 5645                                                                                                                |
| Immunogen        | Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 1-100 of human PRSS2 (P07478). |
| Source           | Rabbit                                                                                                              |
| Category         | Monoclonal Antibodies                                                                                               |
| Application      | WB, IHC-P, ELISA                                                                                                    |
| Cross Reactivity | Human, Mouse, Rat                                                                                                   |
| SWISS            | P07478                                                                                                              |
| Protein Weight   | 26kDa                                                                                                               |
| Shipping         | Ice bag                                                                                                             |