

Anti-Prss23 Polyclonal Antibody (Rabbit Pab) For Wb, Ihc-P, Elisa - O95084

Item Number: CM0022774

Introduction

High-quality rabbit polyclonal antibody against human PRSS23 (serine protease 23), validated for WB, IHC-P, ELISA. Cross-reacts with mouse and rat. Ideal for protease research and IVD assay development.

[Learn More](#)

Parameter	Value
Product Name	PRSS23 Rabbit pAb
Remarks/Alias	SIG13; SPUVE; ZSIG13; PRSS23
Species	Human
Gene ID	11098
Immunogen	Recombinant fusion protein containing amino acids 101-220 of human PRSS23 (NP_009104.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O95084
Protein Weight	43kDa
Shipping	Ice bag