


Anti-Prss8 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q16651

Item Number: CM0020737



Introduction

Rabbit polyclonal antibody against human PRSS8 (prostasin/CAP1), validated for WB, IF/ICC, ELISA. Reacts with human, mouse, rat. Suitable for studying ENaC regulation and epithelial biology.

[Learn More](#)

Parameter	Value
Product Name	PRSS8 Rabbit pAb
Remarks/Alias	CAP1; PROSTASIN; PRSS8
Species	Human
GeneID (Human)	5652
GeneID	5652
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 200 to the C-terminus of human PRSS8 (NP_002764.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q16651
Protein Weight	36kDa
Shipping	Ice bag