

Anti-Psd95 Rabbit Monoclonal Antibody For Wb, Elisa - P78352

Item Number: CM0002904



Introduction

Rabbit monoclonal antibody against PSD95, validated for Western blot and ELISA. Recognizes the N-terminal region of human PSD95 and cross-reacts with mouse and rat. Suitable for detecting endogenous PSD95 in brain tissue, enabling studies of synaptic plasticity and postsynaptic organization.

[Learn More](#)

Parameter	Value
Product Name	PSD95 Rabbit mAb
Remarks/Alias	MRD62; PSD95; SAP90; SAP-90
Species	Human
Gene ID (Human)	1742
Gene ID	1742
Immunogen	A synthetic peptide corresponding to amino acids 1-100 of human PSD95 (P78352).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P78352
Protein Weight	80 kDa
Shipping	Ice bag