

Anti-Psma1 Rabbit Monoclonal Antibody For Wb, If/Icc, Elisa - P25786

Item Number: CM0012144



Introduction

Rabbit monoclonal antibody against human PSMA1, suitable for WB, IF/ICC, and ELISA. Detects human, mouse, and rat PSMA1. Ideal for studying proteasome function and protein degradation pathways.

[Learn More](#)

Parameter	Specification
Product Name	PSMA1 Rabbit mAb
Remarks/Alias	NU; HC2; PROS30; HEL-S-275; PSMA1
Species	Human
GeneID (human)	5682
GeneID	5682
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 164-263 of human PSMA1 (P25786).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P25786
Protein Weight	30kDa
Shipping	Ice bag