


Anti-Psma2 Polyclonal Antibody For Wb, Elisa - P25787

Item Number: CM0017979



Introduction

Rabbit polyclonal antibody targeting PSMA2 (Proteasome subunit alpha type-2), validated for WB and ELISA. Reacts with human and mouse. 26kDa. Ideal for proteasome and ubiquitin-proteasome pathway research.

[Learn More](#)

Parameter	Value
Product Name	PSMA2 Rabbit pAb
Remarks/Alias	MU; HC3; PSC2; PMSA2; PSMA2
Species	Human
Gene ID (Human)	5683
Gene ID	5683
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-234 of human PSMA2 (NP_002778.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P25787
Protein Weight	26kDa
Shipping	Ice bag