


Anti-Psma4 Rabbit Polyclonal Antibody For Wb, Elisa - P25789

Item Number: CM0016613



Introduction

High-quality rabbit polyclonal antibody targeting human PSMA4 (Proteasome subunit alpha type-4). Validated for WB and ELISA. Cross-reacts with human, mouse, rat, and monkey. Ideal for proteasome and protein degradation research.

[Learn More](#)

Parameter	Value
Product Name	PSMA4 Rabbit pAb
Remarks/Alias	HC9; PSC9; HsT17706; PSMA4
Species	Human
GeneID (Human)	5685
GeneID	5685
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-261 of human PSMA4 (NP_001096137.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat, Monkey
SWISS	P25789
Protein Weight	29kDa
Shipping	Ice bag