

Anti-Psma/Folh1 Polyclonal Antibody For Wb And Elisa - Q04609

Item Number: CM0020447



Introduction

Rabbit polyclonal antibody targeting human PSMA/FOLH1 (FOLH1 gene), validated for Western blot (WB) and ELISA. Cross-reacts with human, mouse, and rat. Suitable for prostate cancer biomarker research and glutamate signaling studies.

[Learn More](#)

Parameter	Value
Product Name	PSMA/FOLH1 Rabbit pAb
Remarks/Alias	PSM; FGCP; FOLH; GCP2; PSMA; mGCP; GCP11; NAALAD1; PSMA/FOLH1
Species	Human
GeneID (Human)	2346
GeneID	2346
Immunogen	Synthetic peptide
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q04609
Protein Weight	84kDa
Shipping	Ice bag