

Anti-Psma/Folh1 Rabbit Polymab® For Wb, Ihc-P, Elisa - Q04609

Item Number: CM0011808



Introduction

PSMA/FOLH1 Rabbit PolymAb® is a monoclonal antibody for WB, IHC-P, and ELISA. Detects human, mouse, and rat PSMA, a prostate cancer biomarker and glutamate carboxypeptidase, validated with synthetic peptide immunogen.

[Learn More](#)

Parameter	Value
Product Name	PSMA/FOLH1 Rabbit PolymAb®
Remarks/Alias	PSM; FGCP; FOLH; GCP2; PSMA; mGCP; GCPII; NAALAD1
Species	Human
GeneID (Human)	2346
GeneID	2346
Immunogen	Synthetic peptide corresponding to amino acids 50-150 of human PSMA/FOLH1 (NP_004467.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q04609
Protein Weight	84kDa
Shipping	Ice bag