


Anti-Ptpn13 Polyclonal Antibody For Wb, If/Icc, Elisa - Q12923

Item Number: CM0016272



Introduction

Rabbit polyclonal antibody targeting human PTPN13 (FAP-1/PTPL1), a tyrosine phosphatase involved in apoptosis and PI3K signaling. Validated for WB, IF/ICC, ELISA with cross-reactivity to mouse and rat.

[Learn More](#)

Parameter	Value
Product Name	PTPN13 Rabbit pAb
Remarks/Alias	PNP1; FAP-1; PTP1E; PTPL1; PTPLE; PTP-BL; hPTP1E; PTP-BAS; PTPN13
Species	Human
GeneID (Human)	5783
GeneID	5783
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 250-500 of human PTPN13 (NP_542416.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q12923
Protein Weight	277kDa
Shipping	Ice bag