


Anti-Ptpn4 (Ptpmeg) Rabbit Polyclonal Antibody For Wb, Elisa - P29074

Item Number: CM0028618



Introduction

High-quality rabbit polyclonal antibody targeting human PTPN4 (PTPMEG), suitable for WB and ELISA. Cross-reacts with human, mouse, and rat. Ideal for studying phosphatase signaling, immunity, and neuronal homeostasis.

[Learn More](#)

Parameter	Value
Product Name	PTPN4 Rabbit pAb
Remarks/Alias	MEG; PTPMEG; PTPMEG1; PTPN4
Species	Human
GeneID (human)	5775
GeneID	5775
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-240 of human PTPN4 (NP_002821.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P29074
Protein Weight	106kDa
Shipping	Ice bag