


Anti-Ptrf Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q6Nzi2

Item Number: CM0019007



Introduction

PTRF (Cavin-1) rabbit polyclonal antibody validated for WB, IHC-P, IF/ICC, ELISA. Reacts with human, mouse, rat. Ideal for caveolae formation, adipocyte biology, and metabolic disease research.

[Learn More](#)

Parameter	Value
Product Name	PTRF Rabbit pAb
Remarks/Alias	CGL4, PTRF, CAVIN, FKSG13, cavin-1
Species	Human
Gene ID (Human)	284119
Gene ID	284119
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 200-280 of human PTRF (NP_036364.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q6Nzi2
Protein Weight	43kDa
Shipping	Ice bag