

Anti-Rab11Fip1 Polyclonal Antibody For Wb, Elisa - Q6Wkz4

Item Number: CM0027453



Introduction

Rabbit polyclonal antibody targeting human RAB11FIP1, validated for WB and ELISA. Ideal for studying endosomal recycling, phagocytosis, and vesicle trafficking in research applications.

[Learn More](#)

Parameter	Value
Product Name	RAB11FIP1 Rabbit pAb
Remarks/Alias	RCP; NOEL1A; rab11-FIP1; RAB11FIP1
Species	Human
Gene ID (Human)	80223
Gene ID	80223
Immunogen	Recombinant protein: Recombinant fusion protein containing a sequence corresponding to amino acids 330-540 of human RAB11FIP1 (NP_079427.4).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q6WKZ4
Protein Weight	137kDa
Shipping	Ice bag