

Anti-Rab12 Rabbit Monoclonal Antibody [Ko Validated] For Wb And Elisa - Q6Iq22

Item Number: CM0004941



Introduction

KO validated anti-Rab12 rabbit monoclonal antibody (clone MJFF-A-6E6) for Western blot and ELISA. Detects human and mouse Rab12, a small GTPase regulating endosomal trafficking and autophagy. Ideal for membrane trafficking studies.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] Rab12 Rabbit mAb [MJFF-A-6E6]
Remarks/Alias	RAB12
Species	Human
GeneID (Human)	201475
Immunogen	Recombinant protein containing a sequence corresponding to amino acids 1-244 of human Rab12 protein (Q6IQ22).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q6IQ22
Protein Weight	27kDa
Shipping	Ice bag