

Anti-Rab17 Rabbit Polyclonal Antibody For Wb, Elisa - Q9H0T7

Item Number: CM0026393

Introduction

High-quality rabbit polyclonal antibody against human RAB17, suitable for WB and ELISA applications. Detects RAB17 in human, mouse, and rat samples. Ideal for transcytosis and membrane trafficking studies.

[Learn More](#)

Parameter	Value
Product Name	RAB17 Rabbit pAb
Remarks/Alias	RAB17
Species	Human
GeneID (Human)	64284
GeneID	64284
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-212 of human RAB17 (NP_071894.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9H0T7
Protein Weight	23kDa
Shipping	Ice bag