


Anti-Rab26 Rabbit Polyclonal Antibody For Wb, Elisa - Q9ULw5

Item Number: CM0024448



Introduction

High-quality rabbit polyclonal antibody against human RAB26, validated for WB and ELISA. Reacts with human, mouse, and rat. Targets a key regulator of vesicular trafficking, ideal for neuroscience and cell biology research.

[Learn More](#)

Parameter	Value
Product Name	RAB26 Rabbit pAb
Remarks/Alias	V46133; RAB26
Species	Human
GeneID (Human)	25837
GeneID	25837
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-60 of human RAB26 (NP_055168.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9ULW5
Protein Weight	28kDa
Shipping	Ice bag