


Anti-Rab38 Polyclonal Antibody For Wb, Elisa - P57729

Item Number: CM0022288



Introduction

Polyclonal antibody raised against human RAB38 (NY-MEL-1) for Western Blot and ELISA. Cross-reacts with human, mouse, and rat; detects a 24 kDa band. Useful in melanosome biogenesis and intracellular trafficking studies.

[Learn More](#)

Parameter	Value
Product Name	RAB38 Rabbit pAb
Remarks/Alias	rrGTPbp; NY-MEL-1; RAB38
Species	Human
GeneID (Human)	23682
GeneID	23682
Immunogen	Recombinant Protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P57729
Protein Weight	24kDa
Shipping	Ice bag