


Anti-Rac3 Rabbit Polyclonal Antibody For Wb, Elisa - P60763

Item Number: CM0026539



Introduction

RAC3 Rabbit Polyclonal Antibody (PAb) for WB and ELISA applications. Validated in human and mouse samples. Targeting the small GTPase RAC3 involved in cell spreading and actin dynamics.

[Learn More](#)

Parameter	Value
Product Name	RAC3 Rabbit pAb
Remarks/Alias	RAC3
Species	Human
Gene ID (Human)	5881
Gene ID	5881
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-192 of human RAC3 (NP_005043.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P60763
Protein Weight	21kDa
Shipping	Ice bag