

Anti-Raidd/Cradd Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - P78560

Item Number: CM0004846



Introduction

Recombinant rabbit monoclonal antibody targeting human RAIDD/CRADD (P78560). Validated for WB, IHC-P, and ELISA; cross-reacts with rat. Ideal for apoptosis and cell death research.

[Learn More](#)

Parameter	Value
Product Name	RAIDD/CRADD Rabbit mAb
Remarks/Alias	MRT34; RAIDD; RAIDD/CRADD
Species	Human
Gene ID (Human)	8738
Gene ID	8738
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 70-150 of human CRADD (P78560).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human,Rat
SWISS	P78560
Protein Weight	23 kDa
Shipping	Ice bag