

Anti-Ralb Polyclonal Antibody For Wb, Elisa - P11234

Item Number: CM0030698



Introduction

High-quality rabbit polyclonal antibody targeting human RALB (Ras-related protein Ral-B). Validated for Western blot and ELISA; cross-reacts with human, mouse, and rat. Ideal for studying GTPase signaling, exocytosis, apoptosis, and cell migration.

[Learn More](#)

Parameter	Specification
Product Name	RALB Rabbit pAb
Remarks/Alias	RALB
Species	Human
GeneID (Human)	5899
GeneID	5899
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-206 of human RALB (NP_002872.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P11234
Protein Weight	23kDa
Shipping	Ice bag