


Anti-Rar α Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P10276

Item Number: CM0014517



Introduction

Rabbit anti-RAR α polyclonal antibody validated for WB, IF/ICC, ELISA. Reacts with human, mouse, rat. Targets retinoic acid receptor alpha (P10276), a nuclear receptor involved in transcription, development, and inflammation. Suitable for signaling and epigenetic research.

[Learn More](#)

Parameter	Value
Product Name	RAR α Rabbit pAb
Remarks/Alias	RAR; NR1B1; RARalpha; RAR α
Species	Human
Gene ID (Human)	5914
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-457 of human RAR α (NP_001019980.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P10276
Protein Weight	51kDa
Shipping	Ice bag