


Anti-Rb Rabbit Mab For Wb, Ip, Ihc-P, Elisa - P06400

Item Number: CM0008503



Introduction

High-quality rabbit monoclonal antibody targeting human Retinoblastoma-associated protein (RB1/pRb) for WB, IP, IHC-P, ELISA. Ideal for cell cycle, cancer, and chromatin research.

[Learn More](#)

Parameter	Value
Product Name	RB Rabbit mAb
Remarks/Alias	RB; pRb; OSRC; pp110; p105-Rb; PPP1R130; p110-RB1
Species	Human
Gene ID (Human)	5925
Gene ID	5925
Immunogen	Recombinant protein (Recombinant fusion protein containing a sequence corresponding to amino acids 701-928 of human RB - NP_000312.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IP, IHC-P, ELISA
Cross Reactivity	Human
SWISS	P06400
Protein Weight	106kDa
Shipping	Ice bag