

Anti-Rbbp4 Polyclonal Antibody For Wb, Elisa, Chip - Q09028

Item Number: CM0016578



Introduction

Validated rabbit polyclonal antibody targeting human RBBP4 (NURF55/RBAP48), a histone-binding protein involved in chromatin remodeling. Suitable for WB, ELISA, ChIP. Ideal for epigenetics research.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] RBBP4 Rabbit pAb
Remarks/Alias	NURF55; RBAP48; lin-53; P4
Species	Human
GeneID	5928
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-425 of human RBBP4 (NP_005601.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA, ChIP
Cross Reactivity	Human
SWISS	Q09028
Protein Weight	48kDa
Shipping	Ice bag