

Anti-Rbm19 Rabbit Polyclonal Antibody For Wb, Elisa - Q9Y4C8

Item Number: CM0019475



Introduction

High-quality rabbit polyclonal antibody targeting human RBM19, validated for WB and ELISA. Essential for studying embryo development, RNA-binding mechanisms, and intestinal neoplasia.

[Learn More](#)

Parameter	Specification
Product Name	RBM19 Rabbit pAb
Remarks/Alias	Mrd1; RBM19
Species	Human
GeneID	9904 (Human)
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 759-960 of human RBM19 (NP_057280.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q9Y4C8
Protein Weight	107kDa
Shipping	Ice bag