

Anti-Rbm34 Rabbit Polyclonal Antibody For Wb, Elisa - P42696

Item Number: CM0013884

Introduction

RBM34 rabbit polyclonal antibody validated for Western blotting and ELISA. Targets human RNA-binding protein 34, which localizes to the nucleolus and is implicated in RNA processing. Suitable for research use.

[Learn More](#)

Parameter	Value
Product Name	RBM34 Rabbit pAb
Remarks/Alias	RBM34
Species	Human
Gene ID	23029
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 60-280 of human RBM34 (NP_055829.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P42696
Protein Weight	49kDa
Shipping	Ice bag