


Anti-Recoverin Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - P35243

Item Number: CM0010559



Introduction

High-quality anti-Recoverin rabbit monoclonal antibody validated for WB, IHC-P, ELISA. Suitable for human, mouse, and rat samples. Recoverin is a calcium sensor in phototransduction. Ideal for vision research.

[Learn More](#)

Parameter	Value
Product Name	Recoverin Rabbit mAb
Remarks/Alias	RCV1; RCVRN
Species	Human
Gene ID (Human)	5957
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-200 of human Recoverin (NP_002894.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P35243
Protein Weight	23 kDa
Shipping	Ice bag