

Anti-Relb Rabbit Monoclonal Antibody For Wb, Ihc-P, Ip, Elisa - Q01201

Item Number: CM0007374



Introduction

The RELB Rabbit mAb is a recombinant monoclonal antibody validated for Western blot, IHC-P, immunoprecipitation, and ELISA. It recognizes human, mouse, and rat RELB (62 kDa), enabling studies of the NF- κ B transcription factor pathway.

[Learn More](#)

Parameter	Specification
Product Name	RELB Rabbit mAb
Remarks/Alias	IREL; I-REL; IMD53; REL-B; RELB
Species	Human
Gene ID (Human)	5971
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 296-436 of human RELB (NP_006500.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q01201
Protein Weight	62kDa
Shipping	Ice bag