

Anti-Relb Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q01201

Item Number: CM0016674



Introduction

Rabbit polyclonal antibody targeting human RELB (Transcription factor RelB), validated for WB, IHC-P, ELISA. Cross-reacts with mouse and rat. Ideal for studying NF-kappa-B signaling, lymphocyte development, and inflammation. MW ~62 kDa.

[Learn More](#)

Parameter	Value
Product Name	RELB Rabbit pAb
Remarks/Alias	IREL; I-REL; IMD53; REL-B; RELB
Species	Human
GeneID (human)	5971
GeneID	5971
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 296-436 of human RELB (NP_006500.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB,IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q01201
Protein Weight	62kDa
Shipping	Ice bag