

Anti-Reptin/Ruvbl2 Polyclonal Antibody For Wb, Elisa - Q9Y230

Item Number: CM0013230



Introduction

Rabbit polyclonal antibody targeting Reptin/RUVBL2 (RUVBL2), validated for WB and ELISA in human, mouse, and rat. Involved in chromatin remodeling, transcription, and DNA repair; suitable for NuA4 and INO80 complex studies.

[Learn More](#)

Parameter	Value
Product Name	Reptin/RUVBL2 Rabbit pAb
Remarks/Alias	RVB2; TIH2; ECP51; TIP48; CGI-46; ECP-51; INO80J; REPTIN; TIP49B; TAP54-beta; Reptin/RUVBL2
Species	Human
GeneID (Human)	10856
GeneID	10856
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-463 of human Reptin/RUVBL2 (NP_006657.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y230
Protein Weight	51kDa
Shipping	Ice bag