


Anti-Rgap1 Monoclonal Antibody For Wb, Elisa - Q9H0H5

Item Number: CM0003317



Introduction

RGAP1 Rabbit monoclonal antibody validated for Western blot (WB) and ELISA. Recognizes human RGAP1 (Rac GTPase-activating protein 1), a key cytokinesis regulator. Suitable for cell cycle, cancer, and hematopoiesis research.

[Learn More](#)

Parameter	Value
Product Name	RGAP1 Rabbit mAb
Remarks/Alias	CYK4; CDAN3B; ID-GAP; HsCYK-4; MgcRacGAP; RGAP1
Species	Human
Gene ID (Human)	29127
Gene ID	29127
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-100 of human RGAP1 (Q9H0H5)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q9H0H5
Protein Weight	71 kDa
Shipping	Ice bag