


Anti-Rhoa Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P61586

Item Number: CM0013576



Introduction

CamelBio's anti-RhoA rabbit polyclonal antibody (P61586) is validated for WB, IHC-P, IF/ICC, and ELISA. Detects human, mouse, rat RhoA (22 kDa), a small GTPase regulating cytoskeleton and cell migration.

[Learn More](#)

Parameter	Value
Product Name	RhoA Rabbit pAb
Remarks/Alias	ARHA; ARH12; RHO12; EDFAOB; RHOH12; RhoA
Species	Human
GeneID (Human)	387
GeneID	387
Immunogen	Synthetic peptide
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P61586
Protein Weight	22 kDa
Shipping	Ice bag