


Anti-Rhoa/Rhob/Rhoc Monoclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q16695 / P68431

Item Number: CM0012061



Introduction

Rabbit monoclonal antibody against RHOA/RHOB/RHOC, validated for WB, IHC-P, IF/ICC, ELISA. Cross-reacts with human, mouse, rat. 22kDa protein – Q16695 / P68431.

[Learn More](#)

Parameter	Value
Product Name	RHOA/RHOB/RHOC Rabbit mAb
Remarks/Alias	ARH9; ARHC; H9; RHOH9; RHOA/RHOB/RHOC
Species	Human
Gene ID (Human)	387/388/389
Gene ID	389/387/388
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-193 of human RHOA/RHOB/RHOC (P08134)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q16695/P68431
Protein Weight	22kDa
Shipping	Ice bag