

Anti-Rhod Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - O00212

Item Number: CM0028433



Introduction

Rabbit polyclonal antibody against human RHOD (RhoD), validated for WB, IHC-P, IF/ICC, ELISA. Cross-reacts with mouse and rat. Predicted molecular weight 23 kDa.

[Learn More](#)

Parameter	Value
Product Name	RHOD Rabbit pAb
Remarks/Alias	Rho; ARHD; RHOM; RHOHP1; RHOD
Species	Human
GeneID (human)	29984
GeneID	29984
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 101-200 of human RHOD (NP_055393.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O00212
Protein Weight	23kDa
Shipping	Ice bag