

Anti-Rhog Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P84095

Item Number: CM0023061



Introduction

Rabbit polyclonal anti-RHOG antibody validated for Western blot (WB), IHC-P, IF/ICC, and ELISA. Cross-reacts with human, mouse, and rat RHOG. Ideal for immunological synapse, cell migration, and cancer research.

[Learn More](#)

Parameter	Value
Product Name	RHOG Rabbit pAb
Remarks/Alias	ARHG; RHOG
Species	Human
GeneID (Human)	391
GeneID	391
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 50-150 of human RHOG (NP_001656.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P84095
Protein Weight	21kDa
Shipping	Ice bag