


Anti-Ring1B/Rnf2 Monoclonal Antibody For Wb, Ihc-P, Elisa - Q99496

Item Number: CM0007007



Introduction

Rabbit monoclonal antibody against RING1B/RNF2 (RING2), validated for WB, IHC-P, ELISA. Detects human, mouse, and rat RING1B. Ideal for epigenetic, chromatin, and cancer research.

[Learn More](#)

Parameter	Value
Product Name	RING1B/RNF2 Rabbit mAb
Remarks/Alias	BAP1; DING; BAP-1; HIP13; RING2; LUSYAM; RING1B; RING1B/RNF2
Species	Human
Gene ID (Human)	6045
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-336 of human RING1B/RNF2 (NP_009143.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q99496
Protein Weight	38kDa
Shipping	Ice bag