

Anti-Rip3 Rabbit Monoclonal Antibody For Wb, If/Icc, Elisa - Q9QZL0

Item Number: CM0005582



Introduction

Anti-RIP3 rabbit monoclonal antibody for mouse necroptosis and apoptosis research. Validated for western blot, immunofluorescence/immunocytochemistry, and ELISA applications. Recombinant immunogen from amino acids 320-420 of mouse RIP3.

[Learn More](#)

Parameter	Value
Product Name	RIP3 Rabbit mAb
Remarks/Alias	Rip3; 2610528K09Rik
Species	Mouse
GeneID (Human)	11035
GeneID	56532
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 320-420 of mouse RIP3 (NP_064339.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Mouse
SWISS	Q9QZL0
Protein Weight	63kDa
Shipping	Ice bag